

**Learn more about ReSound Unite™  
wireless accessories.**

**Ask your local hearing specialist  
Scan here or go to  
[www.resound.com/unite](http://www.resound.com/unite)**



[resound.com/  
unite](http://resound.com/unite)



[gnresound.co.uk/  
unite](http://gnresound.co.uk/unite)



[gnresound.com.au/  
unite](http://gnresound.com.au/unite)



[gnresound.co.nz/  
unite](http://gnresound.co.nz/unite)

ReSound® helps people rediscover hearing. We develop hearing solutions that combine original thinking and design with solid technology – enabling you to live a more rich, active and fulfilling life.

ReSound Unite™. A family of proven wireless accessories that extend the reach of your hearing instrument, and complete your hearing system.

## ReSound

red·iscover hearing

The trademarks listed are owned and used by The GN ReSound Group and its related affiliates. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google, Inc.


Let life inspire you  
with ReSound Unite™ wireless accessories

M101528-GB-12.10-Rev.A



## ReSound

red·iscover hearing



Just hit the road

"I have all the accessories,  
and I love them!"

To help you enjoy true hearing freedom, ReSound Unite™ wireless accessories make it easier to engage and interact in any hearing situation.

Our accessories stream sound directly to your hearing aids, so you don't need to wear bulky necklaces. All you need are a few handy devices – and a zest for life.



### ReSound Unite™ Mini Microphone

The meal was great, and for the first time the conversation was just as good. You were even able to continue your chat in the car.

This is no longer wishful thinking. With the ReSound Unite™ Mini Microphone, you don't have to worry about the noise around you when you're out dining, or think about where sounds are coming from. The ReSound Unite Mini Microphone streams voices and sounds directly to your hearing aids – so you can just relax, and enjoy.



## Hear speech clearly

“...changed our life. We can converse in a car and he participates so much more because he can really hear a lot better. We'd never realized how much it had isolated him socially.”

“This is THE BEST technology improvement in many years... In the past, I had to stop what I was doing so that I could read lips. Now I can hear and understand without looking up from the task at hand.”

## Stream audio directly

“I love to listen to baseball games on my...transistor radio! I plug the Mini Microphone in the ear phone plug... If they hit a homerun, I'm grinning, and my wife doesn't know why.”

“I plugged it into our iPod yesterday and just left it while I worked nearby in the yard. I made it at least 30 feet [10 meters] before the signal was affected.”



### ReSound Unite™ Mini Microphone

Streams sound directly to your hearing aids. Elegant and discreet, you simply place it in front of the sound source or clip it on to the speaker you wish to hear. Even plug it into an iPod touch to take your music with you.

- Streams sound at least 7 meters with clear line of sight
- 8½ hours of talk time per battery charge
- 3-hour charge time
- Private and secure connection
- Has microphone for voice pickup and plug for external audio sources

A woman with short brown hair, wearing a dark jacket and a thick, light-colored knitted scarf, stands on a train platform. She is looking down at her smartphone. In the background, a blue high-speed train is stopped at the platform. Other people are visible walking on the platform. A clock and a sign with the number '1' are visible on the wall.

Go mobile

Even when you're in a busy, crowded place, talking on the phone becomes effortless with the ReSound Unite™ Phone Clip+. Just push the answer button and the Phone Clip+ will stream the call directly into both hearing aids.

If it's too noisy, you can even block out surrounding noise to better hear your phone call\*. When you aren't on the phone, use the Phone Clip+ as a basic hearing aid remote control.



### ReSound Unite Phone Clip+

No larger than a matchbox, the Phone Clip+ makes phone calls easier anywhere. Plus it's very handy and easy to use. Clip it on your clothes and keep on talking.

\*Mute function available only with ReSound Verso hearing instruments





“Now I have quality phone calls. I wear the Phone Clip whenever I have my cell phone on.”



[www.resound.com/control](http://www.resound.com/control)



## ReSound Control™

Get your app for iPhone and other smartphones. It gives an easy overview and remote control of your hearing aids and wireless streamers



[resound.com/control](http://resound.com/control)



[gnresound.co.uk/control](http://gnresound.co.uk/control)



[gnresound.com.au/control](http://gnresound.com.au/control)



[gnresound.co.nz/control](http://gnresound.co.nz/control)




### With the ReSound Unite™ Phone Clip+, you get


- 6 hours of talk time and 80 hours of stand-by time
- Call reject and transfer, last number dialed and voice dial functions
- A mute button to block out surrounding sounds so you can concentrate on your call\*
- Simultaneous connection to two Bluetooth® devices (such as phone and iPad)
- Easy-to-use remote control functionality for your ReSound wireless hearing aids
- Quality stereo sound when streaming music\*
- Private and secure connection
- Access to ReSound Control™, an app for iPhone and other smartphones to discreetly control your hearing aid.

Available only with ReSound Verso™ hearing instruments

“Technology that improves my quality of life is paramount”



Discreet  
control  
at your  
fingertips



With the ReSound Unite™ Remote Control, you can easily control all your hearing devices. Its visual interface gives you a clear overview of the current volume and program settings of your hearing aids.

Overview, control and discretion – all in one.



### ReSound Unite™ Remote Control

“For me, it would be awkward to reach behind my ear to control the volume... When I'm in a meeting, I put it on the [correct program] to control the noise.”

#### With the ReSound Unite Remote Control, you can

- Adjust volume of hearing aids together or individually
- Change hearing aid programs
- Start streaming programs and control volume
- Use the lock key to avoid accidental button activation
- Use the home key to take you back to the start screen



Sit back  
and enjoy

Before, you turned up the volume to hear the TV or stereo, and your family complained and left the room. Now, let's invite them back – to a room with comfortable volume and you as a fully engaged conversation partner.

With the ReSound Unite™ TV streamer, you can enjoy clear sound from the TV or any other audio device – at normal volume – and still participate in conversations around you.



### ReSound Unite™ TV

It connects directly to virtually any audio device – TVs, stereos or computers – and delivers entertainment directly to your hearing aids. It's very easy to use and you won't even notice it in your room.



# A world of opportunities

ReSound wireless hearing aids and ReSound Unite™ accessories give you an amazing hearing system so you can get the best out of life.

“My favorite is the TV, because I’ve always had to use subtitles. Now, I can hear everything said, and sung... I can actually make dinner in the kitchen and hear the TV!”

“I had a very hard time understanding TV, even with a hearing aid. And everyone complained that the TV was too loud. Now, my wife can mute the TV, and I can still hear it.”



TV



Mini Microphone



Mini Microphone



Phone Clip+



TV



Phone Clip+



Mini Microphone



Remote Control



Mini Microphone



Mini Microphone



Phone Clip+



Control app



## ReSound Unite™ TV

- Streams stereo sound to both hearing aids
- Streams sound at least 7 meters
- Can be set so you can hear both audio devices and sounds around you
- Can pair and share with as many hearing aids as desired
- Install and forget - always stays paired
- Can easily take it with you to enjoy TV and stereo anywhere